



409 examination table is an economical solution for general medical examinations and minor operations. Thanks to its clear structure 409 is easy and ergonomic to use and quick to clean. Wide range of accessories guarantees comfort and efficiency for medical professionals.

- electric and hydraulic models available
- optional Trendelenburg
- easy-to-clean upholstery
- several colour options: blue, beige, green, bluish grey, grey, black
- push bars at both ends

Accessories and Options

10000732	paper roll holder 732
10000734	paper cutter 734
10000583	wire basket 583, lower frame (hydraulic model only)
100005841	ABS utility basket 5841, head end of the lower frame (hydraulic model only)
100005851	ABS utility basket 5851, leg end of the lower frame (hydraulic model only)
100005861	ABS utility basket 5861 to lower frame for el. models
10000700	accessory rail 700, pair
10000977	Trendelenburg and anti-Trendelenburg adjustment

Technical specifications

Mattress base	2-sectional, epoxy powder coated steel frame
Weight	ca. 80 kg
SWL	180 kg
Length	2030 mm
Width	705 mm
Height	590 - 960 mm
Castors	125 mm, individual braking system 125 mm and 150 mm, central braking system

Adjustments

Backsection	-20 ... +54°, electric or gas spring assisted
Height	590 - 960 mm electric or hydraulic adjustment
Trendelenburg	0 ... 14°, anti-Trendelenburg 0 ... 7°

Warranty 3 years



CE

Provided with specified accessories conforms to the European Council Directive 93/42/EEC, Medical Device Directive Annex VII and the corresponding Finnish National Law No. 1505 (1994) and regulation no. 1506 (1994).

The products is controlled in accordance with ISO 9001 and ISO 13485 Quality Management Systems and ISO 14001 Environmental Management System.